

UPCERA <i>...where art meets science</i>	HT (High Translucent)	ST (Super Translucent)			TT-LT(TT one) (Top Translucent -Low)			TT-MT (Top Translucent - Medium)			Explore Multi-Layered		
	Raw Material	3Y-TZP (3YSB-E)	3Y-TZP (Zpex)			4Y-TZP (Zpex 4)			5Y-TZP (Zpex Smile)			3Y-TZP to 4Y-TZP	4Y-TZP to 5Y-TZP
Restorations	Framework	Full Contour / Monolithic			Full Contour / Monolithic			Full Contour / Monolithic (Aesthetic)			Full Contour / Monolithic		Full Contour / Monolithic (Aesthetic)
Shade	White	White	Pre-shaded (VITA Classical Shades - 16)	Multi-Layered	White	Pre-shaded (VITA Classical Shades - 16/Bleach	Multi-Layered	White	Pre-shaded (VITA Classical Shades - 16)	Multi-Layered	White	Multi-Layered	Multi-Layered
3 Points Bending Strength	1400MPa	> 1200MPa	> 1100MPa	> 1100MPa	> 1000MPa	> 1000MPa	> 1000MPa	> 600MPa	> 600MPa	> 600MPa	800MPa-1100MPa (Incisal - Cervical)		900MPa-600MPa (Incisal - Cervical)
Visible Light Transmittancy	39%	43%	-	-	47%	-	-	49%	-	-	47%-43% (Incisal - Cervical)		49%-46% (Incisal - Cervical)
Indications	Single Crown, Connected Crowns, Bridges (Anterior ≤3 pontics; Posterior ≤2 pontics; 1 pontic in free end bridge)	Single Crown, Connected Crowns, Bridges (Anterior ≤3 pontics; Posterior ≤2 pontics; 1 pontic in free end bridge)			Single Crown, Connected Crowns, Bridges (Anterior & Premolar ≤2 pontics; Posterior ≤1 pontics; 1 pontic in free end bridge)			Single Crown, Connected Crowns, Anterior & Premolar bridges up to 4 units; 3 units Posterior bridges			Single Crown, Connected Crowns, Bridges (Anterior ≤3 pontics; Posterior ≤2 pontics; 1 pontic in free end bridge)		Single Crown, Connected Crowns, Bridges (Anterior ≤1 pontics; Posterior ≤2 pontics; 1 pontic in free end bridge)